Intended use: This chamber is intended to be used by patients who are prescribed a metered dose inhaler or soft mist inhaler, to deliver the aerosol medication from the inhaler to your lungs.

The clinical benefit of chambers is to improve the coordination of inhaling your medication. This is achieved by holding the medication, allowing you more time to inhale while filtering out the larger particles that would not reach your lungs.

If you have any concerns regarding the use of the chamber or its delivery of your inhaler medication, please consult a healthcare professional.

BEFORE EACH USE

This anti-static chamber can be used directly out-of-package. Ensure these instructions and the instructions supplied with the inhaler have been read and are kept available at all times. Carefully look into the chamber. If there is any visible dust or other debris inside, then clean the chamber. Replace immediately if it is damaged or has missing parts.

FOR USE WITH METERED DOSE INHALERS (PUFFERS)

For preparation of the inhaler, follow the instructions supplied with the inhaler.

1 Remove caps from the inhaler and chamber. Shake the inhaler immediately before use as per the instructions supplied with it. **2 MASK:** Insert the inhaler into the backpiece of the chamber. Apply mask to face and ensure an effective seal. **MOUTHPIECE:** Insert the inhaler into the backpiece of the chamber. Put mouthpiece into mouth and close your lips around it to ensure an effective seal. **3 MASK:** Exhale then press the inhaler once at the beginning of a slow inhalation. Hold mask in place and breathe in and out through the chamber for 5-6 breaths. **MOUTHPIECE:** Exhale then press the inhaler once at the beginning of a slow inhalation. Inhale slowly and deeply through the chamber until a full breath has been taken. Hold your breath for 5-10 seconds before exhaling. **OR** Exhale and press the inhaler once at the beginning of a slow inhalation. Breathe in and out through the chamber for 2-3 breaths keeping lips sealed around chamber mouthpiece.

Notes

- Inhalation — SLOW DOWN if you hear the whistle sound, it means you are inhaling too quickly (on adult models).
- If someone is helping you they can use the **Flow-Vu*** Inhalation Indicator to ensure a good seal, coordinate pressing the inhaler with inhalation and to count the number of breaths taken. The **Flow-Vu*** Indicator moves toward you as you inhale and only moves if you have a proper seal.
- Administer one (1) puff at a time. Follow instructions supplied with the inhaler on how long to wait before repeating.
- Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State.

Cautions

- Do not use with Dry Powder Inhalers or Nebulizers. This chamber has not been designed to deliver these medications so you may not get the proper dose.
- Failure to follow instructions may affect medication delivery or cause inhalation issues.

Warnings

- Do not leave the device unattended with children. Contains small parts that can cause choking.

FOR USE WITH RESPIMAT® SOFT MIST INHALERS

For preparation of the inhaler, follow the instructions supplied with the inhaler.

1 Hold the inhaler upright and turn the clear base in the direction of the arrows on the label until it clicks. (Keep the inhaler cap closed). **2** Open the cap of the inhaler until it snaps fully open. Remove the cap from the chamber. **3** Insert inhaler into the backpiece of the chamber. **Do not cover the air vents (A)** on the inhaler. **MASK:** Apply mask to face and ensure an effective seal. **MOUTHPIECE:** Place mouthpiece into mouth and close lips around it to ensure an effective seal. **4 MASK:** Breathe out gently and press the **dose release button (B)** on the inhaler at the beginning of a slow inhalation. Maintain seal for 5 – 6 breaths after the inhaler button is pressed. **MOUTHPIECE:** Exhale then press the **dose release button (B)** on the inhaler at the beginning of inhalation. Inhale slowly and deeply through the chamber until a full breath has been taken. Hold your breath for 5-10 seconds before exhaling. **OR** Exhale then press the **dose release button (B)** on the inhaler at the beginning of a slow inhalation. Breathe in and out through the chamber for 2-3 breaths keeping lips sealed around chamber mouthpiece.

CLEANING

This chamber can be used right out-of-package and then cleaned weekly. (*Mouthpiece model shown as an example*)

1 Remove the backpiece. To detach the frontpiece, twist chamber as shown. **2** Soak the parts for 15 minutes in a mild solution of liquid dish detergent and lukewarm clean water. Agitate gently. Rinse parts in clean water. **OR** Place parts in top rack of dishwasher. Ensure product is securely placed face up. Run the dishwasher on a normal or light cycle. Do not heat dry. **3** Shake out excess water and allow to air dry in a vertical position. Ensure parts are dry before reassembly. **4** To reassemble, fit the frontpiece on the end of the chamber and twist firmly until securely locked into position. Center the alignment feature on the backpiece with the **Flow-Vu*** Inhalation Indicator, as shown. Press firmly to attach the backpiece.

Notes

- Product should be replaced after 12 months of use.
- Product is free of Bisphenol A, Phthalates, Latex, Lead and PVC.
- If you notice medication build-up in your chamber, wash the inside of the chamber gently with a soft cloth.
- Dishwashing with overly dirty dishes is not recommended.
- If cleaning in a dishwasher use a rinse aid.
- For cleaning of the inhaler, follow the instructions supplied with the inhaler.
- The chamber device can be disposed of with domestic waste unless this is prohibited by the disposal regulations prevailing in the respective member countries.
- Put cap back on the mouthpiece and store in a clean, dry area when not in use.
- Testing has been conducted with the following disinfectant solutions:
 - Milton⁺ • Dodie⁺ • Perasafe⁺
- Test results indicate that the solutions tested have no adverse effect on the performance of the device.

Caution

- Do not boil or sterilize. Product may be permanently damaged if boiled, sterilized or cleaned in a dishwasher at a temperature above 70°C.

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


EN Join the myAERO® Community

- Content to help you get started and get the most from your device
- Tips to manage your respiratory symptoms
- Opportunities to participate in contests and special offers
- Annual reminders to replace your device

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AeroChamber
Plus **Flow-Vu**

TRUDELL MEDICAL
INTERNATIONAL

Trudell Medical International
725 Baronsway Drive, London, Ontario,
Canada N5V 5G4
☎ +1-519-455-4862 ☎ +1-519-455-6478
info@trudellmed.com www.trudellmed.com

CE UK
CA

REF 108501 REF 108502 REF 108504
REF 108507 REF 108508 REF 108509

AERO
www.aerochamber.com

103133-001 A, 2024-09-20, * trade marks and registered trade marks of Trudell Medical International (TMI). † trade marks and registered trade marks of their respective companies. © TMI 2024. All rights reserved.

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