
AeroEclipse^{*} XL


BAN^{*} Nebulizer


TRUDELL MEDICAL
INTERNATIONAL

	EN	AEROECLIPSE* XL BAN* Nebulizer	2
	4	المُرذاذ AeroEclipse* XL BAN*	AR
	DA	AEROECLIPSE* XL BAN* Forstøver	7
	DE	AEROECLIPSE* XL BAN* Vernebler	10
	EL	AEROECLIPSE* XL BAN* Νεφελοποιητής	13
	ES (ES)	AEROECLIPSE* XL BAN* Nebulizador	17
	ES (LA)	AEROECLIPSE* XL BAN* Nebulizador	20
	ET	AEROECLIPSE* XL BAN* Nebulisaator	23
	FI	AEROECLIPSE* XL BAN* Nebulisaattori	26
	FR (CA)	AEROECLIPSE* XL BAN* Nébuliseur	29
	FR (FR)	AEROECLIPSE* XL BAN* Nébuliseur	32
	35	נבולייזר AeroEclipse* XL BAN*	HE
	IT	AEROECLIPSE* XL BAN* Nebulizzatore	38
	KO	AEROECLIPSE* XL BAN* 네블라이저	42
	LT	AEROECLIPSE* XL BAN* Nebulaizeris	44
	LV	AEROECLIPSE* XL BAN* Smidzinātājs	47
	MS	AEROECLIPSE* XL BAN* Penebula	50
	NL	AEROECLIPSE* XL BAN* Vernevelaar	53
	NO	AEROECLIPSE* XL BAN* Forstøver	56
	PL	AEROECLIPSE* XL BAN* Nebulizator	59
	PT	AEROECLIPSE* XL BAN* Nebulizador	62
	RO	AEROECLIPSE* XL BAN* Nebulizator	66
	SI	AEROECLIPSE* XL BAN* නිහාරිකාරකය	69
	SL	AEROECLIPSE* XL BAN* Nebulizator	72
	SV	AEROECLIPSE* XL BAN* Nebulisator	75
	TA	AEROECLIPSE* XL BAN* நடுபொலசைர	78
	TH	AEROECLIPSE* XL BAN* เครื่องพ่นยา	82
	TR	AEROECLIPSE* XL BAN* Nebulizatör	85
	ZH (HK)	AEROECLIPSE* XL BAN* 霧化器	88
	ZH (TW)	AEROECLIPSE* XL BAN* 噴霧器	90

1

INTENDED USE

The nebulizer is intended to be used to deliver aerosol drug to the patient's lungs. It is to be used by infant, child and adult patients. The patient is to be awake and can breathe unassisted. The environments of use include hospitals, clinics and the home. The nebulizer is single patient, reusable when used in the home. It is multi-patient, reusable when cleaned according to the hospital instructions.

Breath actuation only produces aerosol when the patient inhales. Aerosol is not produced when the patient exhales or takes a break in treatment. The clinical benefit is more drug is available to the lungs. Less drug is lost to the environment.

⚠ Caution: Breath actuation delivers drug on inspiration only. This may result in an increase of the dosage to the lungs. Health care providers should take this into consideration with regards to dosage.

If you have concerns regarding the use of the nebulizer, please speak with your health care provider.

i Notes

- The nebulizer is designed to aerosolize solutions and suspensions.
- Only use with drugs prescribed by a health care provider.
- Do not use in an anaesthetic or a ventilator breathing system.

BEFORE EACH USE

Ensure these instructions have been read and are kept available at all times. Carefully examine the device. If there is any visible dust or other debris inside, then clean it. Check inside the **EZ Twist Tubing (F)** end for fraying or debris, clean it out if needed. Replace the nebulizer and the tubing **(F)** immediately if they are damaged or have missing parts.

⚠ Warning: Do not leave the nebulizer unattended with children. Contains small parts that can cause choking and tubing that can cause strangulation.

INSTRUCTIONS FOR USE

1 Ensure the exhalation valve **(B)** on the mouthpiece **(A)** faces down.

⚠ Warning: It is important to prevent exhaled breath from blowing near the eyes. Refer to the drug product monograph for more information.

i Note: For patients unable to use the mouthpiece, the nebulizer can be used with a **ComfortSeal[®]** Mask.

2 Inhale and exhale through the nebulizer mouthpiece **(A)** to verify the movement of the green feedback button **(C)** on the top of the nebulizer. The button should move up and down freely.

i Notes:

- To operate the nebulizer in continuous mode **(I)**, rotate the mode selector **(G)** clockwise to the continuous mode position **△**.
- To return the nebulizer to breath actuated mode **(H)**, rotate the mode selector **(G)** counterclockwise to the breath actuated position **⋮**.

3 Unscrew and remove the top of the nebulizer. Place prescribed drug into the nebulizer cup **(E)**. Reattach the top and gently hand tighten.

i Note: Keep the nebulizer in an upright position so the drug does not spill.

4 Use With Compressor: Attach one end of the tubing **(F)** to the fitting in the bottom of the nebulizer and the other end to the compressor. Ensure both ends are securely attached. Turn on the compressor.

i Note: The nebulizer may be used with an **Ombra[®]** Compressor or an air source capable of delivering 103 to 170 kPa (15 to 25 psi) and 4 L/min when connected to the nebulizer.

5 Place the mouthpiece **(A)** in mouth, and inhale slowly and deeply. As you inhale, the green feedback button **(C)** will move into the down position. This indicates the nebulizer is producing aerosol. Exhale normally. As you exhale, the green feedback button **(C)** will move into the up position. This indicates the nebulizer is not producing aerosol.

5 Place the mouthpiece **(A)** in mouth, and inhale slowly and deeply. As you inhale, the green feedback button **(C)** will move into the down position. This indicates the nebulizer is producing aerosol. Exhale normally. As you exhale, the green feedback button **(C)** will move into the up position. This indicates the nebulizer is not producing aerosol.

i Notes

- Do not place the lips over the exhalation valve **(B)** on the bottom of the mouthpiece **(A)**. Doing so will prevent the valve from functioning properly.
 - Ensure the nebulizer remains upright for proper nebulization.
- 6** Continue to inhale from the nebulizer until it starts to sputter, or as instructed by your health care provider.

7 Turn off the compressed air source.

i Note: Disconnect the tubing **(F)** from the compressed air supply after use.

HOME USE – CLEANING

⚠ Caution: Follow the cleaning instructions. Failure to do so may affect drug delivery.

Nebulizer: The nebulizer (excluding the tubing **(F)**) can be cleaned after each use.

1 The mode selector **(G)** should be in the breath actuated position **⋮**.

2 Remove the tubing **(F)** from the bottom of the nebulizer.

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4 Shake out excess drug.

5 Use any of the following methods to clean the nebulizer:

A Hand Washing

i. Soak the three parts in a mild solution of liquid dish detergent and lukewarm clean water for about 5 minutes.

ii. Agitate gently.

iii. Rinse thoroughly with clean water. Shake out excess water.

B Dishwasher

i. Secure the three parts in a basket on the top rack.

ii. Run the dishwasher on a normal cycle. Do not hand dry.

iii. Shake out excess water.

i Note: Dishwashing with overly dirty dishes is not recommended.

6 Briefly connect the tubing **(F)** to the bottom of the nebulizer cup **(E)** and the compressed air supply to clear the air pathway.

7 Allow the parts to air dry thoroughly before reassembling.

8 Reattach the nebulizer top (**D**), making sure to line up the white cylinder with the stem of the nebulizer cup (**E**). Gently hand tighten until a stop is detected.

9 Reattach the mouthpiece (**A**) to the nebulizer with the exhalation valve (**B**) facing down.

10 After the nebulizer is fully dry and assembled, store in a clean plastic bag or container.

EZ Twist Tubing: To remove condensation on the inside of the tubing (**F**), connect one end of the tubing (**F**) to the compressed air supply. Allow the air to run without the nebulizer attached. This will force air through the tubing, drying the inside. To clean the exterior, wipe with a damp cloth. Allow to air dry thoroughly.

HOME USE – DISINFECTION

The nebulizer (excluding the tubing (**F**)) can be disinfected each day.

Follow the cleaning instructions. Do not dry or reassemble the nebulizer.

1 Use any of the following methods to disinfect the nebulizer.

A Boiling

i. Place the three parts in a pot with boiling water and boil for 5 to 15 minutes.

i Note: Place a small rack on the bottom of the pot to prevent parts from coming in contact with the hot bottom.

ii. Remove from the water without touching the water.

iii. Allow the parts to air dry thoroughly before reassembling.

B Electronic Steam Disinfector

i. Place the three parts in an electronic steam disinfector. Follow the manufacturer's instructions for use.

ii. Remove the parts from the unit. Allow the parts to air dry thoroughly before reassembling.

C Microwave Steam Bag

i. Place the three parts in a microwave steam bag (such as Quick Clean[®] Micro-Steam[®] bag). Follow the manufacturer's instructions for use.

ii. Remove from the bag. Allow the parts to air dry thoroughly before reassembling.

D Isopropyl Alcohol

i. Soak the three parts in 70 % isopropyl alcohol for 5 minutes.

ii. Rinse the parts thoroughly with sterile water.

iii. Allow the parts to air dry thoroughly before reassembling.

E Bleach

i. Soak the three parts in a 1:50 bleach to water solution (1 tablespoon (15 mL) bleach in 3/4 cups (800 mL) water) for 3 minutes.

ii. Rinse the parts thoroughly with sterile water.

iii. Allow the parts to air dry thoroughly before reassembling.

Note: You may substitute a respiratory equipment disinfectant (such as Milton[®], Dodie[®], Control III[®] or Cidex[®]) for the bleach solution. Follow the disinfectant instructions for use.

F Vinegar

i. Soak the three parts in a 1:3 vinegar to water solution (1 cup (250 mL) vinegar in 3 cups (750 mL) water) for 1 hour.

ii. Rinse the parts thoroughly with sterile water.

iii. Allow the parts to air dry thoroughly before reassembling.

2 After the nebulizer is fully dry and assembled, store in a clean plastic bag or container.

i Note: The nebulizer (excluding tubing) can be cleaned and disinfected up to 728 times.

HOSPITAL INSTRUCTIONS FOR MULTI-PATIENT RE-USE

⚠ Warning: Clean, disinfect and sterilize the nebulizer before using it with a different patient. This will prevent cross infection.

i Notes:

- The nebulizer (excluding the tubing (**F**)) can be cleaned, disinfected and sterilized up to 150 times.
- The tubing must be replaced between patients.

CLEANING AND DISINFECTION

1 The mode selector (**G**) should be in the breath actuated position .

2 Remove the tubing (**F**) from the bottom of the nebulizer.

3 Unscrew and remove the nebulizer top (**D**).

4 Shake out excess drug.

5 Soak the three parts in a 0.5 % solution of neodisher[®] MediClean forte (or similar) cleaning agent for 5 minutes, agitating gently.

6 Load the parts into a washer-disinfector compliant to ISO 15883-1. Ensure the parts are oriented upright.

7 Run the cleaning and disinfection cycle.

i Notes:

- Use cleaning agent neodisher[®] MediClean.
- The use of a rinse agent is optional.

8 Following disinfection, remove the parts from the washer-disinfector. Shake out excess water. Allow the parts to air dry thoroughly.

STERILIZATION

1 After cleaning, disinfection, and air dry, pack and seal the nebulizer parts in a sterilization pouch compliant to ISO 11607-1 & 2 or EN 868-5.

i Notes:

- The pouch should be no smaller than 70 mm x 120 mm.
- The white cylinder portion of the nebulizer top (**D**) should be positioned in the pouch so it is not bent or pressed into corners where it can become bent (Figure **1**).

2 Load the pouch in the sterilizer.

3 Sterilize according to ISO 17665-1 Moist Heat with any of the following parameter sets:

A

- 3 pre-vacuum phases
- 134 °C sterilization temperature
- 3 minutes hold time
- 1 minute dry time

B

- 3 pre-vacuum phases
- 121 °C sterilization temperature
- 20 minutes hold time
- 1 minute dry time

i Note: Do not exceed a temperature of 137 °C.

4 Reassemble the nebulizer prior to next use.

5 Store in a clean plastic bag or container.

Notes

- No modifications of the nebulizer are permitted.
- The nebulizer parts comply with the relevant requirements of ISO 27427.
- The materials used for the components may not be compatible with emulsions. In such cases, information should be sought from the drug supplier.
- Not made or manufactured with bisphenol A (BPA), di(2-ethylhexyl) phthalate (DEHP), latex or lead.
- This nebulizer can be disposed of with domestic waste unless this is prohibited by the disposal regulations prevailing in the respective member countries.
- Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State, the MHRA or the Saudi Food and Drug Authority, as appropriate.

Type	Pneumatically powered jet nebulizer
Breath Actuation Operation	Triggered mechanically (by the patient on inhalation) or by an inhalation flow over 12 L/min
Driving Gas	Air
Driving Gas Pressure – Compressor (at nominal 4 L/min flow when connected to the nebulizer)	Maximum: 172 kPa (25 psi) Nominal: 140 kPa (20 psi) Minimum: 103 kPa (15 psi)
Driving Gas Flow – Wall Air (at nominal 345 kPa (50 psi) pressure when connected to the nebulizer)	Maximum/Nominal: 8 L/min Minimum: 7 L/min
Maximum Pressure for Tubing at 22 °C ± 2 °C Ambient Conditions	345 kPa (50 psi)
Maximum Fill Volume	6 mL
Maximum A-Weighted Sound Pressure Level	63 dBA
Mass of Nebulizer (without tubing)	0.03 kg (0.07 lb)
Operating Temperature Limits	+15 °C (59 °F) to +40 °C (104 °F)
Operating Relative Humidity Range	15 % to 95 % RH
Storage Temperature Range	-40 °C (-40 °F) to +60 °C (140 °F)
Storage Relative Humidity Range	Ambient to 95 % RH
Recommended Replacement Period	6 months

COMPRESSOR — ISO 27427 PERFORMANCE DATA	At 3.5 L/min (100 kPa)	At 4 L/min (130 kPa)	At 4.5 L/min (160 kPa)
Aerosol Output ¹	10.1 mg	10.8 mg	11.5 mg
Aerosol Output Rate ¹	1.7 mg/min	1.9 mg/min	2.2 mg/min
Percentage of Fill Volume Emitted Per Minute ¹	3.5 %	3.9 %	4.4 %
Residual Volume ¹	1.2 mL	1.1 mL	1 mL
Mass Median Aerodynamic Diameter (MMAD)	4.0 µm ²		
Geometric Standard Deviation (GSD)	2.1 ²		
Particle size % > 5 µm	39 % ²		
Particle size % 2 µm to 5 µm	41 % ²		
Particle size % < 2 µm	20 % ²		
Respirable Fraction (5 µm and smaller)	61 % ²		

WALL AIR — ISO 27427 PERFORMANCE DATA	At 7 L/min	At 8 L/min
Aerosol Output ¹	14.8 mg	16.1 mg
Aerosol Output Rate ¹	3.4 mg/min	3.9 mg/min
Percentage of Fill Volume Emitted Per Minute ¹	6.9 %	7.8 %
Residual Volume ¹	1 mL	1 mL
Mass Median Aerodynamic Diameter (MMAD)	3.4 µm ³	
Geometric Standard Deviation (GSD)	2.3 ³	
Particle size % > 5 µm	33 % ³	
Particle size % 2 µm to 5 µm	40 % ³	
Particle size % < 2 µm	27 % ³	
Respirable Fraction (5 µm and smaller)	67 % ³	

¹ 2 mL fill volume of sodium fluoride 2.5 % (M/V) solution. TMI data on file.

² 2.5 mL fill volume of salbutamol 0.1 % (M/V) solution. TMI data on file.

³ 3 mL fill volume of sodium fluoride 2.5 % (M/V) solution. TMI data on file.

Notes

- A suspension and/or high viscosity solution can alter the particle size distribution curve, the mass median aerodynamic diameter (MMAD), aerosol output and/or aerosol output rate which then can be different from those disclosed in these instructions.
- Disclosures for nebulizer performance are based upon testing that utilizes adult ventilatory patterns and are likely to be different from those stated for paediatric or infant populations.

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