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## EN

### INTENDED USE

The filter kit is intended for use as a single patient, single use accessory to the **AeroEclipse® BAN®** Nebulizer. The clinical benefit of filters is to reduce exposure to potentially harmful nebulizer environmental emissions. This is achieved by filtering aerosolized medication from the exhaled breath.

Aerosolized medication emissions when assembled with **AeroEclipse® BAN®** Nebulizer are ≤ 1% (TMI Data on File).

① **Note:** Ensure these instructions have been read and are always kept available.

### BEFORE EACH USE

Carefully examine the filter kit removing any foreign objects and replace immediately if it is damaged or has missing parts.

### INSTRUCTIONS FOR USE

1 Replace the mouthpiece of the nebulizer with the assembled filter kit.

① **Note:** Set aside the nebulizer mouthpiece. It can be reinstalled after the filter kit has been removed.
**2** Instruct the patient to tightly seal their lips around the mouthpiece, and inhale and exhale normally through the nebulizer. Verify movement of the *green feedback button* (A) on the nebulizer.

**3** Ensure the filter kit’s *exhalation flow path* (B) is not blocked.

4 Continue treatment using the instructions supplied with the nebulizer.

⚠ **Warning:** Single use only. Replace kit after each use. Risk of infection due to multi-patient use.

### ① Notes

- Do not clean or disinfect the filter kit.
- The filter kit can be disposed of with domestic waste unless this is prohibited by the disposal regulations prevailing in the respective member countries.
- Not made or manufactured with bisphenol A (BPA), phthalates, latex or lead.
- Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State, the MHRA or the Saudi Food and Drug Authority, as appropriate.

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