

# Effect of Inhalation Patterns on the Delivered Dose of Symbicort<sup>†</sup> from a Dry Powder Inhaler Compared to a Metered Dose Inhaler plus Valved Holding Chamber

Am J Respir Crit Care Med 2020;201:A4763  
Authors: M. Nagel, R. Ali, J. Suggett, C. Doyle;  
TruMed Medical International, London, ON, Canada.


## Rationale:

- The latest guidelines from the Global Initiative for Asthma (GINA) no longer recommends SABA-only treatment for Step 1 therapy. This is based on evidence that a SABA-only treatment increases the risk of severe exacerbations, and that adding any Inhaled Corticosteroid (ICS) significantly reduces that risk.
  - Options for ICS delivery can be via Dry Powder Inhaler (DPI) or Metered Dose Inhaler (MDI), however, if delivering an ICS via MDI, a spacer should also be prescribed.
  - This study evaluated the impact of inhalation maneuvers on medication delivery from either DPI or MDI Symbicort<sup>†</sup> 80/4.5 (80 µg budesonide /4.5 µg formoterol furoate) with a spacer (**AeroChamber Plus\* Flow Vu\* VHC**).
- 



# Methods:

- Inhalation waveforms were recorded from 5 DPI and MDI+VHC subjects and replayed via a breathing simulator attached to the adult Aerosol Delivery to an Anatomic Model (ADAM) oropharyngeal airway which was used to provide a clinically relevant laboratory determination of medication delivery of Symbicort<sup>†</sup>.
- DPI subjects inhaled as per their usual treatment regimen whereas VHC subjects were instructed to inhale via either i) slow deep inhalation or ii) following 2-3 inhalations as per VHC instructions.
- Following delivery of medication, all components of the apparatus were assayed for Active Pharmaceutical Ingredients (API) via HPLC-spectrophotometry.






# Results:


Product	API	Delivered Dose ( $\mu\text{g}/\text{actuation}$ )	
		Oropharyngeal Airway	Carina
Symbicort DPI	Budesonide	21.1 $\pm$ 8.1	20.2 $\pm$ 7.9
	FF	1.0 $\pm$ 0.4	1.0 $\pm$ 0.4
Symbicort MDI + VHC Deep Inhalation	Budesonide	7.1 $\pm$ 3.9	35.9 $\pm$ 4.0
	FF	0.3 $\pm$ 0.2	2.1 $\pm$ 0.2
Symbicort MDI + VHC Tidal Breathing	Budesonide	9.7 $\pm$ 9.6	32.0 $\pm$ 6.8
	FF	0.4 $\pm$ 0.4	1.9 $\pm$ 0.4

## Percentage of Delivered Dose of Budesonide via each Delivery Method





## Conclusion



- The widely differing inspiratory flow profiles from Turbuhaler<sup>†</sup> DPI contributed to greater variability in delivered dose to the carina.
- MDI delivery via the **AeroChamber Plus\* Flow-Vu\*** VHC resulted in significantly increased mass to the carinal region irrespective of inhalation maneuver as well as reduced oropharyngeal deposited mass which *in vivo* is likely to result in reduced throat irritation or thrush in the mouth and throat.