

Complete Delivery of Medication from a Valved Holding Chamber is Affected by Chamber Capacity

M. Nagel, J. Suggett, V. Avvakoumova
Trudell Medical International, London, Ontario, Canada.

RATIONALE

- Valved Holding Chambers (VHCs) are often prescribed to optimize fine particle delivery of inhaled corticosteroid medication, to the lungs.
- Large volume VHCs have capacities greater than can be emptied in a single breath.
- This laboratory based study compared a widely prescribed large and small volume VHC, for the delivery of Pressurized Metered Dose Inhaler (pMDI)-delivered Fluticasone Propionate (FP).

MATERIALS AND METHODS




Volumatic⁺ VHC with mouthpiece
Large volume 750 mL



AeroChamber Plus* Flow-Vu* Antistatic VHC with mouthpiece
Smaller volume 149 mL

- Each VHC was connected to a breathing simulator (ASL 5000, Pittsburgh, PA) mimicking a tidally breathing adult
 - Tidal volume = 500 mL
 - Rate/min = 13
 - Inspiratory:Expiratory ratio; 1:2
- Aerosol filter collected the medication delivered at the exit of the VHC




- 1 actuation of FP (125 µg) was delivered to the VHC
- Mass collected on the filter after 1 complete breathing cycle ($n=5$ replicates) was assayed by HPLC
- The procedure was repeated after 2, 3, 5, 8 and 10 breathing cycles

- $n=3$ devices/group; 3 replicates/device
- VHCs were prepared per instructions for use
 - Volumatic⁺ — wash in mild detergent, rinse, drip-dry
 - **AeroChamber Plus* Flow-Vu*** VHC — no preparation required

RESULTS

Table 1: Total Mass of pMDI-Delivered FP Recovered

Breathing Cycles	mean \pm SD	
	Volumatic ⁺ VHC	AeroChamber Plus* Flow-Vu* VHC
1	21.9 \pm 2.9	51.2 \pm 4.5
2	14.9 \pm 4.4	53.5 \pm 4.9
3	26.9 \pm 2.8	58.9 \pm 7.8
5	30.0 \pm 4.2	58.8 \pm 4.7
8	29.9 \pm 5.6	54.1 \pm 7.8
10	34.8 \pm 3.9	52.6 \pm 5.4



DISCUSSION

- Differences in particle transport through the two different sizes of VHC investigated may have been partly responsible for the observed behavior.
- Another possible explanation for the lower and erratic performance with the Volumatic⁺ VHC relates to the likely presence of electrostatic charges acquired on the internal surface of that VHC through normal handling even though pre-washing was undertaken.
- Other potential reasons for variability in data may relate to inhalation valve performance.

CONCLUSIONS

- This study highlights that not all chambers perform the same and may deliver different amounts of medication.
- Clinicians should be aware that the volume of a VHC has a marked effect on the delivery of medication.
- An effective VHC should be able to maintain the availability of medication for inhalation in cases where multiple inhalations (and therefore more time) may be required to empty.

