

AeroTrach Plus* Anti-Static Valved Holding Chamber

AeroTrach Plus* Anti-Static Valved Holding Chamber (VHC) is specifically designed for delivery of metered dose inhaler (MDI) medications to tracheostomy patients in hospital and home environments. It can be easily and comfortably connected to standard 15mm tracheostomy tubes.

Refer to the [instructions for use](#) for additional information.

- Use of a chamber improves medication delivery¹
- Anti-static improves the reliability of medication delivery²
- Ready to use directly out of the package – time saving with no pretreatment needed
- Reusable for 4 months

Comfort Right Fit* Adaptor
Attaches and detaches without excessive torque for patient comfort

Transparent Chamber
Clear chamber allows for visualization of aerosol

Inhalation Valve
Easy actuation with minimal inspiratory effort

Exhalation Valve
Low resistance valve minimizes any aerosol loss to the atmosphere between inhalations

Universal Backpiece
Can be used with most commonly prescribed metered dose inhalers



Part Number	Case Quantity
101505110010	10

SINGLE-PATIENT – MULTIPLE USE for 4 months

PACKAGING: Resealable storage bag with instructions for use

CLEANING: Liquid dish detergent and lukewarm clean water

DISINFECTION: Bleach solution, VitalOxide[†] or enzyme solvent

Not made or manufactured with bisphenol A (BPA), phthalates, latex or lead

1. [TARGETING DRUGS TO THE AIRWAYS: THE ROLE OF SPACER DEVICES](#), Lavorini F, et al. Expert Opinion on Drug Delivery 2009;6(1):91-102.

2. [USE OF VALVED HOLDING CHAMBERS WITHOUT PRE-CONDITIONING AND THE INFLUENCE OF ANTI-STATIC MATERIALS](#), Suggett J, et al. Journal of Aerosol Medicine and Pulmonary Drug Delivery 2014;28.

AeroTrach Plus*

Anti-Static Valved Holding Chamber



+1-519-455-4862
info@trudellmed.com



Designed and Developed in Canada
Manufactured in Canada with USA and Canadian parts.

trudellmed.com